

Polymerase Chain Reaction Protocol for Mycoplasma Detection

Materials and reagents:

- 0.2 ml PCR tubes + rack
- pipettes
- ice / cooler
- autoclaved ddH₂O
- reagents:
 - reaction buffer (10 x Optimized DyNAzyme React.Buffer (Thermo Fisher))
 - DNA polymerase (DyNAzyme (Thermo Fisher))
 - dNTPs (10 mM (Sigma))
 - primers (12.5 pM, Mycoplasma forward and reverse (IDT))
 - template DNA sample (DNA from GI261, SW480 cells)

Procedure

1. prepare reaction mixture without DNA for 3 samples:

samples: a. SW 480 (Mycoplasma +)
b. GI261 (Mycoplasma –)
c. blank (= sample without DNA)

mixture: [40 µl final volume / sample]

- water (ddH ₂ O):	27 µl
- 10 x reaction buffer:	5 µl
- 10 mM dNTP:	1 µl
- 12.5 pM Primers:	3-3 µl
- DNA polymerase:	1 µl

for 1 sample: 40 µl

2. add 40 µl reaction mixture into the PCR tubes
3. add DNA into the PCR tubes: 50 ng DNA + ddH₂O (final volume: 10 µl)
4. carefully flick the tubes
5. place them into PCR machine and set the program:

cycle 1	95 °C 3 minutes	
cycle 2	step 1. 95 °C 30 sec	} (40x)
	step 2. 50 – 70 °C 30 sec	
	step 3. 72 °C 1 min	
cycle 3	72 °C 1 min	

hold 4 °C or take the tubes immediately and keep them at 4°C.

Primer sequences: reverse: 5'– TTC TTT TCA CCT TTC CCT CAC GGT AC - 3'

forward: 5'– GGT GAA TAC GTT CTC GGG TCT TGT ACA CAC – 3'