

Mitochondrial DNA Deletion by Real-time PCR Protocol

Reagents and Materials:

- **Template:** 500 ng total DNA
- **Primer set:**
 - Genomic housekeeping gene: GAPDH
 - Mitochondrial housekeeping gene: TotMT
 - WT mitochondrial specific: WT
 - Common deletion specific: DEL
- **Master MIX (MM):** SYBR Green Master mix
- **Prepare Master Mix with the appropriate primer:**

Measure the appropriate amounts of primer with MM in an Eppendorf tube and mix well. Centrifuge briefly and keep ice-cold.

Work flow:

➤ **Primer solution:**

- **MM-primer solution1:**

SYBR Green Master mix: 12.5 µl	}	15.5 µl/tube
primer GAPDH: 3 µl		
- **MM-primer solution2:**

SYBR Green Master mix: 12.5 µl	}	15.5 µl/tube
primer TotMT: 3 µl		
- **MM-primer solution3:**

SYBR Green Master mix: 12.5 µl	}	15.5 µl/tube
primer WT: 3 µl		
- **MM-primer solution4:**

SYBR Green Master mix: 12.5 µl	}	15.5 µl/tube
primer DEL: 3 µl		

➤ **Template solution:**

DNA: x µl (=500 ng DNA)	}	9.5 µl/tube
DNase free water: 9.55 µl – x µl		

Mix the MM primer solution (1 or 2 or 3 or 4) with Template solution

Total volume: 25 µl

- Aliquot the solutions into the 0.2 ml tubes.
- Place the tubes into the 36 rotor. Place the rotor into the thermal cycler.
- Start the RotoGene program.
- New run: mt DEL

PCR program: 10 min. 95 °C
 15 s 95 °C }
 30 s 60 °C } 40x
 30 s 72 °C }

Melt:

95 °C 1 min.

65 °C 2 min.

65-95 °C; 2 °C/min

Save. File name: